

# Outbreaks, Investigations and Evidence for Flexible Endoscope Reprocessing – Overcoming the Challenges

---

HONG KONG HOSPITAL AUTHORITY - ADVANCED COURSE FOR INFECTION CONTROL – NOVEMBER 2017

# Learning Objectives

---

- Present overview of recent endoscope related outbreaks
- Discuss factors that have contributed to the use of contaminated flexible endoscopes
- Understand key recommendations in new endoscope reprocessing guidelines and standards
- Review quality control requirements
- Discuss practical aspects and rationale for using various methods to detect damaged or dirty endoscopes, including enhanced visual inspections, biochemical indicators, and microbial cultures

# Janet Prust - Disclosure

**Employee of 3M Health Care**

**Infection Prevention Division**

Association for Advancement of Medical Instrumentation (AAMI)

Positions held:

AAMI Board of Directors – Director representing industry since 2015

Member:

AAMI Finance committee

Sterilization Standards Committee

WG 61: Chemical sterilants hospital practices – co-chair

WG 84: Endoscope reprocessing

WG 40: Steam sterilization hospital practices

WG 13: Washer disinfectors; TAG to ISO TC 198 WG 13

WG 93: Cleaning of reusable devices

Sterilization of endoscopes stakeholders group

Task group – HVAC conditions in OR




# The Outbreaks

---

# Flexible Endoscope Remain Contaminated After Reprocessing – Alfa 2012

Report of contaminated endoscopes well before publicized outbreaks



# Contamination and Outbreaks Occurring – Kovaleva **2013**

Report of outbreaks occurring in published literature before media reporting

Endoscope	Outbreaks	Primary organism	Patient contaminated	Patients infected	Root cause
Upper GI	19	<i>P. aeruginosa</i> , <i>H. pylori</i> , <i>Salmonella</i>	169	56	Cleaning or disinfection gaps
Sigmoidoscopy/ colonoscopy	5	<i>Salmonella</i> , HCV	14	6	Cleaning or disinfection gaps
ERCP	23	<i>P. aeruginosa</i>	152	89	C/D, water bottle, contaminated AER
Bronchoscopy	51	<i>P. aeruginosa</i> , Mtg, Mycobacteria	778	98	C/D, AER, Water
<b>Totals</b>	<b>98</b>		<b>1,113</b>	<b>249</b>	

## **Reported gastrointestinal endoscope reprocessing lapses: The tip of the iceberg.**

Alexandra M. Dirlam-Langlay, Cori L. Ofstead, Natalie J. Mueller, Pritish K. Tosh, Todd H. Baron, Harry P. Wetzler. American Journal of Infection Control 2013 Dec;41(12):1188-94.

- Looked for reprocessing lapses in peer-reviewed literature, gov't reports, state health depts, CDC, FDA, Dept. of Veteran affairs and media reports
- The study was limited to Jan. 2005 – June 2012.
- ***They found that improper endoscope reprocessing is an ongoing and pervasive problem.***
- Over 30,500 people exposed and this is just the “tip of the iceberg”.
- > 99% of these cases were not found in peer-reviewed medical journals

**Reprocessing lapses are rarely reported in medical journals leading to the false conclusion that reprocessing lapses are rare.**

# First reports of duodenoscope CRE related outbreaks

US government action related to outbreaks occurring in published literature before media reporting

- 2012 outbreaks with multi-drug resistant organisms seen
- US CDC published alert January 2014

The screenshot shows the CDC website interface for a Morbidity and Mortality Weekly Report (MMWR) article. The article title is "Notes from the Field: New Delhi Metallo-β-Lactamase-Producing *Escherichia coli* Associated with Endoscopic Retrograde Cholangiopancreatography – Illinois, 2013". The article is dated January 3, 2014, and is volume 62(51), pages 1051-1051. The text discusses the increasing prevalence of carbapenem-resistant *Enterobacteriaceae* (CRE) in medical facilities, specifically focusing on New Delhi metallo-β-lactamase (NDM)-producing *E. coli* associated with ERCP procedures in Illinois. It details a case-control study conducted at hospital A, showing a strong association between ERCP and the acquisition of NDM-producing CRE. The study found that 8 out of 8 ERCP patients (100%) were case-patients, compared to 4 out of 27 (15%) controls. The odds ratio was 78.0, with a 95% confidence interval of 6.0 to >999.99. The article also mentions that after manual cleaning and high-level disinfection of an automated endoscope reprocessor, cultures were obtained from the ERCP endoscope used on five of the case-patients. NDM-producing *E. coli* and KPC-producing *K. pneumoniae* were recovered from the terminal section (the elevator channel) of the device. The *E. coli* isolate was highly related (>95%) to the outbreak strain by PFGE. Retrospective review and direct observation of endoscope reprocessing did not identify lapses in protocol. Previous studies have shown an association between ERCP endoscopes and transmission of multidrug-resistant bacteria; the design of the ERCP endoscopes might pose a particular challenge for cleaning and disinfection (2,3). Finally, the article states that among 91 ERCP patients who were initially notified that they had potential exposure to a culture-positive endoscope, 50 returned for rectal surveillance cultures. NDM-producing *E. coli* were recovered from 23 (46%). An additional 12 patients with NDM-producing CRE have been identified in northeastern Illinois, bringing the total during January–December 2013 to 44. In September 2013, as a result of the investigation, hospital A changed ERCP endoscope reprocessing from automated high-level disinfection to gas sterilization with ethylene oxide; no new cases with exposure to a gas-sterilized ERCP endoscope have been identified.



# The Outbreaks: 2015 - In the news but not new.....

*After a well documented history of outbreaks*

**NEW YORK** **LOS ANGELES**

**BREAKING OVERNIGHT**  
**SUPERBUG OUTBREAK AT UCLA**  
NEARLY 200 PEOPLE POSSIBLY EXPOSED


**SUPERBUG OUTBREAK**

Caused by contaminated medical scopes

UCLA Health Center

**FDA**  
"QUESTIONING MANUFACTURERS  
ABOUT METHODS FOR CLEANING"

# High levels of Persistent Contamination on Patient Ready Endoscopes – Ofstead 2013, 2014, 2016



Citation: Ofstead et al. The effectiveness of reprocessing in accordance with current guidelines. SGNA Conference Poster. 2015.

# Published CRE Outbreaks – A Wake Up Call

Primary Author	Multi-drug resistant organism	Number of patients impacted	Yr. of publication
Carbonne	<i>K. Pneumoniae</i>	13	2010
Smith	<i>E. Coli</i>	3	2015
Marsh	<i>K. Pneumoniae</i>	34	2015
Kim	<i>K. Pneumoniae</i>	15	2016
Epstein	<i>E. Coli</i>	39	2015
Kola	<i>K. Pneumoniae</i>	12	2015
Wendorf	<i>E. Coli</i>	35	2015
Vertaillie	<i>P. Aeruginosa</i>	22	2015
Total		173	

# Additional Cited Outbreaks:

1. Bajolet O, Ciocan D, et.al. **Gastroscopy-associated transmission of extended-spectrum beta-lactamase-producing *Pseudomonas aeruginosa***. J Hosp. Infect 2013 (83)

---

2. Epstein L, Hunter J, et. Al. **New Delhi Metallo  $\beta$ -Lactamase-producing carbapenem-resistant *Escherichia coli* associated with exposure to duodenoscopes**. JAMA 2014. (312:1447-55)
3. Kovaleva J, Degener J., et.al. **Methylobacterium and its role in health care associated infection**. J Clin Microbiol. 2014 (52). 1317-21
4. Wendelboe A, Baumbach J., et.al. **Outbreak of cystoscopy related infections with *Pseudomonas aeruginosa*: New Mexico. 2007**. J Urol 2008. (N180). 588-92.
5. Wendorf K, Kay M., et.al. **Endoscope retrograde cholangiopancreatography-associated Amp C *Escherichia coli* outbreak**. Infect Control Hosp Epidemiol 2015. (#^). 634-42.
6. Guy M, Vanhems P. , et.al. **Outbreak of pulmonary *Pseudomonas aeruginosa* and *Stenotrophomonas matlophilia* infections related to contaminated bronchoscope suction valves, Lyon, France**. Euro Surveil. 2014. 2016:21.

# Additional Evidence Detail

- ❖ Saliou, et al. 2016. **Persistent Contamination on Endoscopes – Micro surveillance post processing with 8 day incubation**
  - ❖ 34% exceeded target of <25 CFU of indicator organism
- ❖ England, et al. 2016. Transmission of MDRO from **gastroscope** – 5 patients with CRE (9 more exposed)
  - ❖ Superbug persisted through 12 reprocessing cycles
- ❖ **Duodenoscope** ERCP / CRE Outbreaks
  - ❖ Illinois 2013 – 156 patients exposed to CRE/ 39% transmission rate
  - ❖ Washington 2014 - 39 cases/ 18 deaths.
  - ❖ Wisconsin 2014 – 3 deaths
  - ❖ California – 2014 - 2 deaths
- ❖ DiasGranadose, et al. 2009. **Bronchoscope** related outbreak.
  - ❖ 19 patient exposed, 12 infected, 2 deaths
- ❖ FDA MDRO Reports:
  - ❖ 2015. Contaminated **cystoscopes**. 4 patients infections
  - ❖ 2015. Contaminated **ureteroscopes**. 7 infections. Damaged endoscopes.

Outbreaks occurring with all types of flexible endoscopes and a variety of organisms


# Flexible GI Endoscope are not the only risk....

## Reprocessed Flexible Bronchoscopes: FDA Safety Communication - Risk of Infection Sept. 7, 2015

### FDA analysis to date has identified two recurrent themes:

- Failure to meticulously follow manufacturer instructions for reprocessing
- Continued use of devices despite integrity, maintenance and mechanical issues.”

**FDA Recommendation:** “Implement a comprehensive reprocessing quality control program. Your reprocessing program should include written procedures for monitoring, training and adherence to the program, and documentation of equipment tests, processes, and quality monitors used during the reprocessing procedure.”



# Other reprocessing lapses and patient exposure to contaminated endoscopes

US Location	Facility types	Errors in reprocessing	Patients impacted
Colorado	Medical center	Improper cleaning	71
Minnesota 1	ASC; Outpatient clinic; 5 hospitals	7 incidents reported: improper cleaning/HLD; reprocessing single use device; inadequate training	6 – 2000 per incident
North Carolina	Hospital	No cleaning/sterilization of one channel	10
New Jersey	ASCs	Improper reprocessing; unchanged water/cleaning solution	Not reported
Ontario, Canada	Clinic	Multiple cleaning/HLD breaches	6800
Louisiana	Medical center	Wrong HLD temperature	360
British Columbia, Canada	Hospital	Bioburden allowed to dry before cleaning	536
California	Hospital; surgery center	Improper HLD; expired disinfectant	3400
Minnesota 2	Medical center	No HLD of one channel	2600
Florida	Hospital; cancer treatment center	Improper cleaning of elevator channel	191
Georgia	Surgery center	Wrong HLD time	1300

# Non-endoscopic Related Outbreaks

---

## Poor Infection Control Practices:

- New York – improper handling of intravenous sedation tubing, multi-dose vials and reuse of needles.
- Las Vegas – Outbreak of Hepatitis C due to cross-contamination from syringes - ~40,000 patient exposed. No follow-up information.



# Outbreaks Related to Tubing, Water Bottles and other sources


---



- Hep C outbreak related to malfunctioning one-way valves contaminated water bottles (water and bottles not changed between patients)
- Biofilm formation on tubing
- Contamination of endoscope from hands of HC workers
- Failure to sterilize biopsy forceps


Ofstead et al., Assessing residual contamination, *AJIC*, 2016;  
Ofstead et al., Longitudinal assessment, *AJIC*, 2017.

**143 gastroscope testing encounters (June-October)**



# The Outbreaks: The microbes are changing the game

- Carbapenem-resistant Enterobacteriaceae – CRE
- Limited or no treatment
- High transmission rate 6-46%
- High mortality rate ~ 50%



THREAT LEVEL  
**URGENT** ○ ○ ○ ○ ○

These bacteria are immediate public health threats that require urgent and aggressive action.

## MICROORGANISMS WITH A THREAT LEVEL OF URGENT

- Clostridium difficile*
- Carbapenem-resistant Enterobacteriaceae
- Drug-resistant *Neisseria gonorrhoeae*



# Guidelines

---



# FDA Safety Communication Feb. 19, 2015

---

- The complex design of duodenoscopes may impede proper reprocessing
- Meticulous manual cleaning should reduce risk of transmission of infection
- Implement a comprehensive Quality Control program
- Quarantine scopes suspected of association with patient infection until shown to be free of pathogens

<http://www.fda.gov/medicaldevices/safety/alertsandnotices/ucm434871.htm>



# Gastroenterology and Urology Devices Panel Meeting May 14-15, 2015

*FDA seeking expert scientific and clinical opinion. 19 member advisory panel.*

“Duodenoscopes and AERs do not provide a reasonable assurance safety and effectiveness”

“Manual Cleaning is a critical component.”

There is a need for “...development and validation of cleaning verification assays”

“Majority of the panel also believes it is necessary to reclassify duodenoscopes from semi-critical to critical and support the move from high level disinfection to sterilization.”

# Supplemental Measures to Enhance Duodenoscope Reprocessing. FDA Safety Communication. August 4, 2015.

<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm454766.htm>

---

## Enhanced measures include:

- Ethylene Oxide Sterilization.
  - Multiple rounds of High-Level Disinfection
  - Use of a liquid chemical sterilant processing system
  - Microbial surveillance
- Implementation of these additional measures **may not be feasible**
- The **limitations** of each of these measures must be taken into consideration.
  - These measures are to be considered **in addition to** following **manufacturer's instructions** for reprocessing, **meticulous manual cleaning** and the implementation of a **comprehensive quality control program**.

This recommendation currently questioned regarding lack of evidence for liquid disinfection recommendations.

# Preventable Tragedies: Superbugs and How Ineffective Monitoring of Medical Device Safety Fails Patients, Jan.13, 2016

United States Senate Health, Education, Labor and Pensions Committee

Patty Murray, Ranking Member

<http://www.help.senate.gov/imo/media/doc/Duodenoscope%20Investigation%20FINAL%20Report.pdf>

Between 2012 and spring 2015, closed-channel duodenoscopes were linked to at least 25 different incidents of antibiotic-resistant infections that sickened at least 250 patients worldwide.

Hospitals, FDA and mfr's *all failed* in their responsibility to report, notify and act on knowledge that outbreaks were occurring.





# New Guidelines Issued

**SGNA (4):** Society for Gastroenterology Nurses and Associates

**AAMI:** Association for Advancement of Medical Instrumentation

**AORN:** Association of periOperative Registered Nurses

**APSIC:** Asia Pacific Society of Infection Control

**ASGE:** American Society for Gastrointestinal Endoscopy/ Multisociety update

**CDC/HICPAC:** Healthcare Infection Control Practices Control Advisory Committee

## To varying level of detail, guidelines agree on:

- Implement quality control program
- Follow manufacturer's instructions for reprocessing
  - New, revalidated IFUs
- Pay special attention to manual cleaning
  - Meticulous cleaning required
  - Elevator mechanism needs special attention
- Implement comprehensive training with competency assessment by qualified personnel
- Periodic review of policies and procedures
- AORN, AAMI, APSIC and SGNA recommend use or assessment of cleaning verification
- ASGE, CDC recommends microbial surveillance

# SGNA Guidelines

- Standards of Infection Prevention in Reprocessing Flexible Gastrointestinal Endoscopes 2015
- Standards in Infection Prevention in Gastroenterology Settings 2015
- Standards of Infection Prevention in Reprocessing Flexible Gastrointestinal Endoscopes 2016
- Guideline for Use of High-Level Disinfectants & Sterilants in the Gastroenterology Setting 2017

<http://www.sgna.org/Portals/0/Standards>

# ANSI/AAMI ST91:2015 Flexible and semi-rigid endoscope reprocessing in health care facilities – **currently under revision**

---

- Design of the endoscope processing area, including work flow considerations
- Personnel issues such as training, hygiene, clothing, policies, and immunizations
- Processing steps and recommendations including sterilization
- Addresses other related topics including automated endoscope reprocessors (AER), sterile endoscope sheaths, and processing accessories
- Storage and transport
- Quality control
- Bibliography



# CDC Interim Microbial Surveillance Protocol -

## March 2015

- Protocol suggested but not yet validated
- Sensitivity of this method is unknown
- Look for pathogens and elevated levels of non-pathogens
- Frequency of testing not defined
  - Weekly, monthly, every time, every 60 procedures
- Pay Special attention to
  - Inspection and Manual Cleaning
  - Drying

<http://www.cdc.gov/hai/outbreaks/index.html>

### Safety Alerts:

#### ERCP Duodenoscopes

- [Interim Duodenoscope Surveillance Protocol](#)
- [Interim Duodenoscope Sampling Method](#)
- [Interim Duodenoscope Culture Method](#)
- [CDC Statement: Los Angeles County/UCLA investigation of CRE transmission and duodenoscopes](#)
- [NDM-Producing CRE Associated with ERCP](#)
- [NDM-Producing CRE Associated With Duodenoscope Exposure](#)
- [See CRE homepage](#)

#### U.S. Food and Drug Administration (FDA) resources:

- [Design of Endoscopic Retrograde Cholangiopancreatography \(ERCP\) Duodenoscopes May Impede Effective Cleaning](#)
- [FDA Guidance: Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling](#)
- [Notice of Meeting: Gastroenterology and Urology Devices Panel of the Medical Devices Advisory Committee](#)

# CDC Interim Protocol: The Jury is Still OUT.....

---

“ ...Not sufficient in the current form to be implemented by healthcare facilities as best practice” *FDA Panel on Gastroenterology and Urology, May 14-15, 2015*

Sensitivity unknown *CDC Interim Protocol for Duodenoscope Surveillance*

“...clinical microbiology labs should not perform routine cultures of reprocessed duodenoscopes due to lack of data on utility of such culturing” *American Society for Microbiology statement on CDC Interim Protocol.*

# Discussion - Outbreaks

---

1. Why are we seeing the increase level of outbreaks?
2. What has changed?
3. How can these be detected?
4. Any experiences in your facilities to share?
5. What is the level of awareness and concern by administration/management? Auditors?



# Human Factors

---

# Ofstead C, Wetzler H, Snyder A. Endoscope Reprocessing Methods: A prospective study on the impact of human factors and automation. Gastroenterology Nursing, 2010. Vol 33 (54): 3-10.

(N = 69 GI endoscopes)

2 or more steps performed improperly nearly half of the time ***while being observed***





# Reprocessing staff comments


---

- ❖ 50% didn't like performing manual cleaning
- ❖ 53% reported discomfort due to working with scopes
- ❖ 13% problems with lungs, sinuses, or breathing
- ❖ 27% bothered by odors or fumes
- ❖ 50% physical discomfort or bothersome symptoms
- ❖ 75% felt pressure to work quickly while reprocessing
- ❖ 37% observed procedural delays due to reprocessing

# IAHCSMM Survey: 2017

Staff report of issues encountered during reprocessing.

N = 67



# Factors that contribute to patient risk related to flexible endoscopy procedures:

---

- ❖ Complicated reprocessing guidelines
- ❖ Lack of training by qualified instructors
- ❖ Non-adherence to guidelines
- ❖ Repeated exposure of endoscopes to organic debris and pathogens
- ❖ Inability to see pathogens and internal endoscope damage with the naked eye
- ❖ Resistance to testing for contamination
- ❖ Increasing pressure to do more, faster, with less resources
- ❖ Longstanding belief that there's no risk – 1 in 1.8 million often cited by GI community

# Flexible Endoscope Reprocessing – Basic Steps\*:

---

- ❖ Bedside pre-cleaning and transport
- ❖  Leak testing
- ❖  Thorough manual cleaning
- ❖  Visual inspection and verification
- ❖  Rinsing
- ❖  High-level disinfection
- ❖  Rinsing
- ❖  Drying
- ❖  Storage



\* This is an oversimplification. It is a very complex process with over 200 steps.

# Complicated Guidelines

Guideline or Standard	Pages	# reference to IFU
AORN Guideline for processing flexible endoscopes (2016)	84	73
SGNA Standards of Infection Prevention in Reprocessing Flexible GI Endoscopes (2015)	31	66
SGNA Standards of Infection Prevention in the Gastroenterology Setting (2015)	22	13
SGNA Guideline for Use of High-Level Disinfectants & Sterilants in Gastroenterology Setting. (2017)	5	13
AAMI ST 91: Flexible and semi-rigid endoscope processing (2015)	70	149
CDC (US Centers for Disease Control and Prevention). Essential Elements of a Reprocessing Program for Flexible Endoscopes – Recommendation of the Healthcare Infection Control Practices Advisory Committee (2017)	12	4
APSIC	9	12

# Clinical Use of products Containing Simethicone

## Active ingredient: silicone ...

- In bowel prep formulations
- Lubricants/de-foaming agents used in procedures, e.g. simethicone
  - Silicone is NOT WATER SOLUABLE
- Add substance to water bottle or injected into endoscopes using a syringe to reduce foam, bubbles, gas
- Sprayed on surfaces to facilitate procedure



- Other ingredients in products include thickeners, sweeteners, preservatives, and water

# Discussion – Human Factors

---

1. Describe your experience regarding staff or human factors issues?
2. What works well? What doesn't.
3. What are the key factors that prevent effective reprocessing from a staff or management factors?
4. What Guidelines are available and followed?
5. How are staff trained and assessed?

# Process Effectiveness and Quality Control

---



# How Fast do Microbes Multiply?

Time	# organisms
After HLD reprocessing	1
20 minutes	2
40 minutes	4
1 hour	8
2 hours	64
3 hours	512
4 hours	4,096
5 hours	32,768
6 hours	262,144
7 hours	2,097,153

Prompt initiation of cleaning is critical to an effective cleaning process.

Current guidelines recommend documentation of time pre-clean was completed to determine if delayed reprocessing procedures are needed.

# Implementation of a Quality Control (QC) Program



*“Implement a comprehensive quality control program for reprocessing duodenoscopes. Your reprocessing program should include written procedures for monitoring training and adherence to the program, and documentation of equipment tests, processes, and quality monitors used during the reprocessing procedure.”*

## **Design of Endoscopic Retrograde Cholangiopancreatography (ERCP) Duodenoscopes May Impede Effective Cleaning: FDA Safety Communication**

What to consider? Quality Control through the implementation of rapid cleanliness indicators, ethylene oxide sterilization to replace HLD and appropriate use of microbial monitoring audits.

# Endoscope Reprocessing Quality System Requirements

---

- ❖ Training and routine competency assessments
  - ❖ Certification
- ❖ **Visual inspection**
- ❖ **Cleaning Verification**
- ❖ Disinfectant concentration testing
- ❖ Physical/mechanical monitoring
- ❖ **Water quality testing to ensure no contamination**
- ❖ Ensure AER / Washer Disinfection functioning and **not contaminated**
- ❖ **Effective drying**
- ❖ **Storage and handling to prevent contamination**
- ❖ **Validation of equipment and processes\***

\* Required in some regions

# Training and Competency

---

- ❖ **Complete and thorough training programs**

- ❖ Standards
- ❖ Manufacturer's instructions for use
- ❖ Certification through formal program recommended

- ❖ **Competency assessment**

- ❖ Direct observation of competency
- ❖ Routine assessment (e.g. quarterly)
- ❖ Assessment performed by knowledgeable, experienced individual in endoscope reprocessing

# Visual Inspection for gross soil and damage

## SGNA

- “Due to recent issues with reprocessing, **visual inspection warrants its own step. It** may be considered a “time out” or safety stop to verify that the endoscope is at least visually clean before proceeding to HLD.”
- “Look for conditions that could impact HLD (corrosion, discoloration, cracks, and retained debris)

## AORN

- “**Lighted magnification should be used** to inspect endoscopes and accessories for cleanliness and damage.”
- “An endoscope that appears clean may harbor debris that cannot be seen without magnification. Lighted magnification may increase the ability to identify residual soil or damage.”
- “Internal channels of flexible endoscopes may be inspected using an endoscopic camera or borescope.”
- “**Defective endoscopes, accessories, and equipment should be removed from service and repaired or replaced**”



Brand new



Damaged




# Borescopes – a new approach to inspection

---



- Small fiberoptic camera to view inside channels
- Greatly magnified image
- Available for S/B channel size
- Can capture images
- Valuable for inspection of damage and gross debris
- Challenges:
  - No benchmark as acceptable wear and tear
  - Disinfection between uses
  - Sizes not available for all channels
- Guidelines referencing for inspection

# Borescope images



# What does everyone agree on?

---

## Focus on Manual Cleaning

- It is a problem
- It is critical to success of HLD or Sterilization
- Lack of proper manual cleaning contributed to outbreaks
- It can be improved
- **Use validated, real-time indicators of cleaning efficacy**
  - Commercially available kits that test for ATP, protein, hemoglobin, carbohydrate





# ANSI/AAMI ST91:2015 Section 12: Quality Control

## QC is critical to successful reprocessing

---

- All facilities should have a comprehensive QC program
  1. Product identification and traceability
  2. Documentation and record-keeping
  3. **Verification and monitoring of the cleaning process**
    - Technologies: ATP, protein, hemoglobin, carbohydrate, bioburden
    - Provides Pass/Fail threshold measuring effective soil removal
  4. Monitoring of high-level disinfection and sterilization processes
  5. Product recalls procedure
  6. Quality process improvement program



# Cleaning Verification Recommendations

## **AORN:**

“Manual cleaning of flexible endoscopes **should be verified using cleaning verification tests** when new endoscopes are purchased and **at established intervals (eg, after each use, daily).**”

## **AAMI:**

“The frequency of **testing the efficacy of the manual cleaning step should occur on a regular basis, weekly or preferably daily...**”

## **SGNA:**

If the test results for organic soil are positive, the endoscope should be **re-cleaned** before HLD

## **APSIC:**

Cleaning verification should be performed at least daily and prefer every time the endoscope is manual cleaned.

# Cleaning Verification Practices


---

## **Samples endoscope surfaces and channels after manual cleaning:**

- Swab the insertion tube, distal end, and ports
- Flush sterile water through channels and collect effluent
- Run rapid biochemical cleaning verification tests for residual contamination
  - ATP
  - Protein (less sensitive)
  - Hemoglobin (only detects presence of blood)
- Re-clean whenever results exceed benchmarks for “clean”
- Re-test after re-cleaning
- Repeat until test passes or a decision is made to send it for repair

# Duodenoscopes Have Added Complexity

- Elevator Guide wire channel: Sealed or Open?
- Elevator housing and mechanism on distal tip.




# Cleaning = Effective Removal of Clinical Soil

## What is in Clinical Soil?


- Microorganisms: bacteria, fungi, viruses
  - Microbes found in all soil components
- Tissue
- Blood and other body fluids
- Secretions/Excretions (vomit, diarrhea, mucous, phlegm, etc)

## Can presence of soil be measured?

- Yes, but it is not easy



# Components Found in Clinical Soil




To measure effectiveness of cleaning - a component is selected to measure.

# Common Universal Cleanliness Markers


Present in all types of soils, used for general monitoring

---

## Protein



## Adenosine Triphosphate (ATP)



adenosine triphosphate (ATP)

# Characteristics of Cleaning Markers

## ATP

Universal Marker

Mature technology

Commercially available

Rapid tests are objective, numeric

ATP is stable under reprocessing conditions

Will detect microbial ATP

Cannot be used when prions are an issue

Different manufacturer's use different measurement scales so cannot compare different systems.

## Protein

Universal Marker

Mature technology

Commercially available

Rapid tests are subjective, colorimetric

Rapid tests cannot measure proteins that are denatured under reprocessing conditions.

Does not detect microbial proteins: requires additional sample prep for detection

Used when prions are an issue

Measurement scale is in standardized  $\mu\text{g}$



# Rapid Cleaning Verification Tests: Rapid Protein Detection

---

- Qualitative Test
- Colorimetric (color change)
- Subjective interpretation
- Gives presence/absence results
- Analysis of data is limited
- Validation data unavailable for most tests
- Data needed to verify sensitivity claims
  - Sensitivity dependent on sample prep?
- Mature technology, well known



# Rapid Cleaning Verification Tests: Multi-test Strip

---

- Three pads with protein, carbohydrate and hemoglobin tests
- Qualitative, subjective interpretation
- Presence/absence test
- Analysis of data is limited
- Validation data shows little utility for protein and carbohydrate pads, hemoglobin test drives detection of dirty scopes
- Because Hg test is specific, the number of dirty scopes detected is lower than expected, when compared to current literature.



# Rapid Cleaning Verification Tests: ATP Bioluminescence

---


- Quantitative numerical results
- Results are given in Relative Light Units (RLU)
- Different mfr's use different measurement scales
- Numerical results can be analyzed – statistics!
- Validation data available (must be carefully interpreted)
- Technology has been used to assess cleaning efficacy for over 35 years.  
(Food Safety, Aerospace, Clean-Room manufacturing)




# ATP Bioluminescence Technology

## Converts ATP to a light signal

Fire-fly enzyme Luciferase uses ATP to produce Light



### Simple Relationship



# Benchmark for Clean

---

- The Pass/Fail benchmarks were developed to assess if a scope has been cleaned (soil removal) to an acceptable level.
  - Acceptable levels for ATP, protein, hemoglobin, carbohydrate, bioburden, endotoxin and sodium ion are defined in AAMI ST91
- *Pass/Fail benchmarks for any cleaning verification test are not a measure of the risk of pathogen transmission, nor does a “pass” guarantee that the scope is safe to use on a patient.*
- Ensure Pass/Fail thresholds have been validated for the technology used.

# Quantitative or Qualitative Cleaning Verification Test

---

## Things to consider

- Is there a validation study available?
- Are there studies that show the product can detect dirty scopes in a clinical setting?
- Do the claims make sense given the technology used?
- Ease of interpretation
- What are you getting for your time and money?



# Standards and Guidelines Positions on Cleaning Verification

---

APUSIC 2017, AAMI ST91 2015, AORN 2016, SGNA 2016 recommend the use/ assess the use of rapid cleanliness indicators for verification of efficacy of manual cleaning.




# Discussion – Quality Control

---

1. Are current quality control practices acceptable for flexible endoscopes?
2. What are the factors that prevent better quality control?
3. Cleaning verification – how often should it be performed?
4. Use of lubrication aids – how does this impact cleaning effectiveness? How can policy be established to reduce/eliminate use?
5. What is needed for better training and competency assessment?



Clinical soil removal



ATP  
Bioluminescence



Microbial  
Surveillance

Effectiveness of entire process

# Microbial Surveillance of Endoscopes

## Guideline recommendations

- ESGE/ ESGENA– Minimum every three months
- ASGE/ Multisociety: Consider monthly for duodenoscopes
- GoSA: Varies based on scope type; monthly for high risk. Included monthly testing of AER and Water
- APSIC: Periodic surveillance

## Periodic surveillance of duodenoscopes recommended by US CDC 2015, 2017

- Controversial
  - Lack of validated method and requires environmental culture
  - Questions on appropriate frequency
  - Questions on validity as method does not detect all pathogens
  - CDC Recommendation is 1x/mo or after 60 procedure

## AORN 2016 and SGNA 2017

- Recommend to conduct a risk assessment and may be considered in event of an outbreak
- Can be used as a method for assessing quality of reprocessing or defects in endoscopes

## AAMI 2015: No guidance; 2018 - TBD

# Specific Tests: Used for Investigations or Answering Very Specific Questions.

## Microbiological Counts

- Viable microbial numbers
- Pathogen detection
- Effectiveness of entire process



## Hemoglobin

- Blood



# Culturing of Endoscopes: Limitations

Current methods are not sensitive enough to detect low-levels of bacteria, limitations of these methods not being discussed

---

## Current methods:

- Do not detect all bacteria
- Do not detect viruses or parasites
- Do not substantiate cleanliness
- Do not substantiate any level of sterilization or disinfection



## Current methods not sufficient for sampling duodenoscopes

- Biofilm bacteria must be cultured differently
- Bacteria exposed to disinfectants need special culture conditions

# Should ATP measurements correlate with bacterial counts?

## *A common concern, a common misconception....*

### Non-culture methods: ATP Bioluminescence


- ATP: measures cleaning effectiveness of reduction in soil to 'threshold' validated value
  - ATP measures organic contamination from all living sources
    - Microorganisms, Human cells, secretions, excretions, body fluids, food residue

### Microbial culture:

- Does not 'correlate' to ATP
- Measures bacteria only
- Indication of effectiveness of entire disinfection or sterilization process

Correlation of ATP and Microbial Counts is not possible as they measure two different things and provide different pieces of Quality Control information.

# Sequence of QC testing



# Recommendations

No longer consider Endoscopy as a low risk procedure


## Implement QA programs

- Written policies
- Training/Competencies
- Regular Audits and continued oversight
- Make sure IFUs are up to date
- Informed consent for patients

## Implement a monitoring program for manual cleaning

- Multiple sampling sites
- Multiple methods
- Every scope, every time

## Implement periodic microbial surveillance for HLD procedures as an audit tool



# Summary Points

---

- Patient-ready endoscopes are contaminated
- Viable microbes commonly survive and patient infection transmission occurs
- Patient outcomes have been catastrophic and public health is at risk due to superbugs
- Guidelines and instructions for use are too complex
- Employees rarely do or are able to follow the steps in guidelines and IFU
- Quality improvement is urgently needed



Thank you for your time and  
attention!

# References - web

- CDC Threat Report: [www.cdc.gov/drugresistance/threat-report-2013/](http://www.cdc.gov/drugresistance/threat-report-2013/)
- FDA Safety Communication. Feb. 19, 2015  
<http://www.fda.gov/medicaldevices/safety/alertsandnotices/ucm434871.htm>
- FDA Supplemental Measures. Aug. 4, 2015  
<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm454766.htm>
- Gastroenterology and Urology Devices Panel Meeting May 14-15, 2015  
<http://www.fda.gov/AdvisoryCommittees/CommitteesMeetingMaterials/MedicalDevices/MedicalDevicesAdvisoryCommittee/Gastroenterology-UrologyDevicesPanel/ucm445590.htm>
- AAMI ST91 Flexible and semi-rigid endoscope reprocessing in health care facilities  
<http://my.aami.org/store/SearchResults.aspx?searchterm=st91&searchoption=ALL>
- CDC Interim Duodenoscope Protocol Mar.12, 2015 <http://www.cdc.gov/hai/outbreaks/index.html>

# References (2)

- **Reprocessed Flexible Bronchoscopes: FDA Safety Communication - Risk of Infection** Sept. 7, 2015  
<http://www.fda.gov/Safety/MedWatch/SafetyInformation/SafetyAlertsforHumanMedicalProducts/ucm462979.htm>
- Alexandra M. Dirlam-Langlay, Cori L. Ofstead, Natalie J. Mueller, Pritish K. Tosh, Todd H. Baron, Harry P. Wetzler. **Reported gastrointestinal endoscope reprocessing lapses: The tip of the iceberg.** American Journal of Infection Control 2013 Dec;41(12):1188-94.
- Ofstead, Cori L., Wetzler, Harry, P., Alycea Snyder, Rebecca A. Horton, **Endoscope Reprocessing Methods: A Prospective Study on the Impact of Human Factors and Automation.** 2010 Gastroenterology Nursing. Vol 33, No. 4, pp. 304-311
- Alfa et al. 2006 American Journal of Infection Control 34, 561-570
- Bommarito, M., Thornhill, G., Morse, D. 2013. **A Multi-site Field Study Evaluating the Effectiveness of Manual Cleaning of Flexible Endoscopes using an ATP Detection System.** American Journal of Infection Control 41(6) p. S24.
- Alfa MJ, Fatima I, Olson N. 2013. **The adenosine triphosphate test is a rapid and reliable audit tool to assess manual cleaning adequacy of flexible endoscope channels.** Am J Infect Control.41:249-53.
- Sciortino CV Jr, Xia EL, Mozee A. 2004. **Assessment of a novel approach to evaluate the outcome of endoscope reprocessing.** Infect Control Hosp Epidemiol. 25:284-90.
- **Bernstein, B. et al. Dramatic Improvement in Endoscope Cleanliness Using a Commercial Bioburden Detection System and Implementation of a Modified Cleaning Algorithm.** Gastrointestinal Endoscopy, Vol.79 (5), Supplement. May 2014, Pages AB121 DDW 2014 ASGE
- Kovaleva, J., Peters, F.T., van der Mei, H.C., & Degener, J.E. (2013). **Transmission of infection by flexible gastrointestinal endoscopy and bronchoscopy.** *Clinical Microbiology Review*, 26(2), 231–254.

# References (3)

- Bajolet O, Ciocan D, et.al. **Gastroscopy-associated transmission of extended-spectrum beta-lactamase-producing *Pseudomonas aeruginosa***. J Hosp. Infect 2013 (83)
- Epstein L, Hunter J, et. Al. **New Delhi Metallo  $\beta$ -Lactamase-producing carbapenem-resistant *Escherichia coli* associated with exposure to duodenoscopes**. JAMA 2014. (312:1447-55
- Kovaleva J, Degener J., et.al. **Methylobacterium and its role in health care associated infection**. J Clin Microbiol. 2014 (52). 1317-21
- Wendelboe A, Baumbach J., et.al. **Outbreak of cystoscopy related infections with *Pseudomonas aeruginosa*: New Mexico. 2007**. J Urol 2008. (N180). 588-92.
- Wendorf K, Kay M., et.al. **Endoscope retrograde cholangiopancreatography-associated AmpC *Escherichia coli* outbreak**. Infect Control Hosp Epidemiol 2015. (#^). 634-42.
- Guy M, Vanhems P. , et.al. **Outbreak of pulmonary *Pseudomonas aeruginosa* and *Stenotrophomonas matlophilia* infections related to contaminated bronchoscope suction valves**, Lyon, France. Euro Surveil. 2014. 2016:21.
- Ofstead CL, Wetzler HP, Heymann OL, Johnson EA, Eiland, JE, Shaw MJ. Longitudinal assessment of reprocessing effectiveness for colonoscopes and gastroscopes: Results of visual inspections, biochemical markers, and microbial cultures. *American Journal of Infection Control*, 2016; In press.
- Ofstead CL, Wetzler HP, Eiland, JE, Heymann OL, Held, SB, Shaw MJ. Assessing residual contamination and damage inside endoscopes over time. *American Journal of Infection Control*, September 2016; Epub ahead of print.
- Ofstead CL, Wetzler HP, Johnson EA, Heymann OL, Maust TJ, Shaw MJ. Simethicone residue remains inside gastrointestinal endoscopes despite reprocessing. *American Journal of Infection Control*, July, 2016; 44(11): 1237-40.
- O'Horo, J. C., Farrell, A., Sohail, M. R., et al. (2016). Carbapenem-resistant Enterobacteriaceae and endoscopy: An evolving threat. *Am J Infect Control* 44, 9: 1032-6.

○